
CDA Lab safety bulletin


November, 2020

Right use of fume hoods

Definition of a Fume Hood



the laboratory worker




Vertical Sash

- ✓ Work with sashes down as low as practical
- ✓ This gives protection against fumes, vapors as well as splashes



Common Sash Types





All fume hoods at KAUST are combination

Horizontal Sash



- ✓ Close the sashes as much as is practical
- ✓ And work with the sash IN FRONT of you
- ✓ If there is a risk of explosion/implosion, use a blast shield in addition to the sash



Golden Rules For CHEMICAL FUME HOODS



- 1. Make sure the hood is working.
- Keep the sash down except when working in the fume hood.



- 3. Raise the sash only the minimum needed to perform the work.
- 4. Work six to eight inches inside the plane of the sash.



- 5. Keep hood area clear and rear vent unobstructed.
- 6. Raise equipment off hood deck. (Use PVC pipe donuts, not wood.)





- 7. Remove flammables and unnecessary stuff.
- 8. Understand the limitations.



- Check / inspect operation at least annually, sash cables every 2-3 years.
- 10. Wear proper personal protective equipment (PPE).






- 11. Run cords, cables and tubing under the air foil.
- 12. Move arms slowly in and out of the hood.

PPE (Personal protective Equipment)

DON'T BE ON THE WRONG SIDE OF SAFETY

Wet Lab Personal Protective Equipment (PPE)



Face coverings are public health measures during the COVID-19 global Pandemic, and not considered personal protective equipment (PPE).

Lab Safety Recognition Awards




If you know someone that meets the following criteria?

- Initiates and/or lead a successful safety initiative;
- Regularly reports unsafe working laboratory conditions and works to correct them;
- Promotes safe behavior and a safe working environment inside the laboratory;
- Leads changes in procedures, training or other administrative controls that help improve safety within the laboratory;
- Engages lab members in creating a safety culture within the lab.

Nominate them for the Laboratory Safety Recognition Award. This is open to anyone working in a laboratory.

Nomination Link: Lab Safety Recognition Awards-2020 - Formstack









SAFETY

SAFETY IS EVERYONE'S RESPONSIBILITY